

Workflow

Batch 1

1. **Add Points:** Left click on objects in frame 0.
2. **Generate Masks:** Press R (Single frame) or click "Single" button
3. **Propagate:** Press M or click "Batch" button to propagate through frames 0-99
4. **Move to Next Batch:** Press N (auto-saves masks)

Batch 2+ (With Previous Masks)

1. **Review Masks:** Previous batch's masks automatically appear on frame 0
At frame 0:
Option A - Good Masks: Press M immediately to propagate
Option B - Need Corrections:
 - a. Add correction points or new points on frame 0.
 - b. Press R to generate new mask for that frame (If new points added)
 - c. Press M to propagate to all frames in batch.
2. Move to a different frame:
 - a. Add new points or corrections.
 - b. Press R to generate new masks for that frame.
 - c. Press M to propagate from current frame to end of batch.
3. **Move to Next Batch:** Press N (auto-saves masks)

Note: If you want to save a batch of masks to a different folder (other than default),

1. Press T.
2. Type the name of the folder in terminal and hit enter.

Controls

Mouse Controls

Action	Description
Left Click	Add positive point (include in mask)
Middle Click	Add negative point (exclude from mask)
Shift + Click	Remove nearest point

Keyboard Shortcuts

Key	Action
A / Z	Switch to previous/next object
; / '	Navigate backward/forward 1 frame
[/]	Jump backward/forward 5 frames
R	Run SAM2 on current frame only (Preview)
M	Merge & propagate from current frame forward (Batch Process)
T	Manually save masks (optional, auto-saves on batch change)
D	Delete last point on current frame
N	Next batch
Q / ESC	Quit application

On-Screen Buttons

- **<<5 / <1 / 1> / 5>>**: Frame navigation
- **Single**: Run SAM2 on current frame (same as R)
- **Batch**: Propagate from current frame forward (same as M)